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【GRUNDINFORMATIONEN】

JOYSBIO(Tianjin) Biotechnology Co., Ltd.

Adresse: Tianjin International Joint Academy of Biotechnology &

Medicine 9th floor No.220, Dongting Road, TEDA 300457 Tianjin

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
EYBEL

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INTENDED USE

For in vitro qualitative detect of SARS-CoV-2 nucleocapsid antigen in nasal(NS) swab specimens directly from individuals who are suspected of COVID-19 by their healthcare provider within the first 5 days of the onset of the symptoms. This test is only provided for use by clinical laboratories or to healthcare workers for point-of-care testing, and not for at home testing.

SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

PRODUCT PHOTOS

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SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

PACKAGE SIZE (BOX)

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Product package size:


Length:195mm

Width:165mm

Height:68mm

Weight:308g

Include:20 Test Kit, 20 pcs/box

SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

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
TEST PRINCIPLE

The Kit use immunocapture method, it is designed to detect the presence or absence of SARS-CoV-2 nucleocapsid proteins in respiratory samples from patients with signs and symptoms of infection who are suspected of COVID-19.


TEST METHOD

1


1. Before collecting oral fluid relax your cheeks and gently massage cheeks with fingers for 15-30 seconds, Gently spit oral fluid into the collection bag.

2

2. Hold the dropper vertically and **draw oral fluid from** collection bag and transfer 3 drops of oral fluid into the buffer bottle.

3


3. Twist off the top of the buffer bottle, slowly dispense all of the buffer into the extraction Tube.


SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

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
4. Tighten the cap of the buffer bottle. Gently shake the buffer bottle for **10 seconds**.

5. Tear off the foil pouch, take out the test strip/cassette and place the test kit on a clean and level surface. Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well. Read the test results between 15 and 20 minutes.

6. **POSITIVE**: Two lines appear. One colored line should be in the control line region (C), a colored line appears in test line (T) region. **NEGATIVE**: Only one colored control line appear. **INVALID**: Control line fails to appear.

CLINICAL EVALUATION REPORT

The kit showed 95.10% of sensitivity and 100% of specificity.

Table 1. Clinical Study Results from symptom onset

Reagent test results	PCR Comparator		Subtotal
	positive	negative	
positive	97	0	97
negative	5	260	265
Subtotal	102	260	362

Positive Percent Agreement (PPA)= $97/102(95.10\%)$
(95%CI:88.9%~98.4%)

Negative Percent Agreement (NPA)= $260/260(100\%)$
(95%CI:98.6%~100%)

Accuracy= $(97+260)/362 \times 100\%=98.62\%$

Kappa= $2 \times 25220 / 52250 = 0.97 > 0.5$

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SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

S I G N I F I C A N C E

SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

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RESEARCH BACKGROUND

During the epidemic Situation, many countries have the following problems:

Existing detection methods cannot achieve large-scale rapid screening.

lack of technical expertise and inadequate laboratory capacity, Erroneous Operation can easily lead to missed inspections.

Can't afford high testing costs.

SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

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SIGNIFICANCE

According to the WHO, during the outbreak of SARS-CoV-2, in areas with confirmed SARS-CoV-2 community-wide transmission; confirmed outbreaks in closed or semi-closed communities; in high-risk groups; among contacts of confirmed cases; SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold) as a tool to monitor disease incidence is a particularly effective detection method.

Globally, as of 3:59pm CEST, 17 August 2020, there have been 21,549,706 confirmed cases of COVID-19, including 767,158 deaths, reported to WHO.

SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

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ADVANTAGE

1. Easy to collect samples, simple operation, without professional equipment.
2. The test results are available in 15 minutes, and the test results are clearly visible.
3. Convenient transportation and low price, higher accuracy.
4. Suitable for large-scale rapid screening.

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Dear Hongyan Li:

The Center for Devices and Radiological Health (CDRH) of the Food and Drug Administration (FDA) has received your submission. This submission has been assigned the unique document control number below. All future correspondence regarding this submission should be identified prominently with the number assigned and should be submitted to the Document Control Center at the above letterhead address. Failure to do so may result in processing delays. If you believe the information identified below is incorrect, please contact the Office of Product Evaluation and Quality (OPEQ) submission support at (301) 796-5640 or OPEQSubmissionSupport@fda.hhs.gov.

Submission Number: EUA202733
Received: 9/11/2020
Applicant: JOYSBIO (Tianjin) Biotechnology Co., Ltd.
Device: SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

We will notify you when the review of this document has been completed or if any additional information is required. For information about CDRH review regulations and policies, please refer to <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/default.htm>.

Sincerely yours,

Center for Devices and Radiological Health

SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

Has entered the FIND recommended list

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FIND  Because diagnosis matters

COVID-19 WHO WE ARE WHAT WE DO NEWSROOM PARTNERS & DONORS

CONTACT US

CALLS FOR PARTNERS

- [Hunan Yonghe-Sun Biotechnology Co., Ltd](#) SARS-CoV-2 specific antibody test kit (Immunochromatography) (RUO) [Contact](#)
- [InDevR Inc.](#) COVID Serology Kit: Multiplexed Immunoassay (RUO) [Contact](#)
- [Innovita Biological Technology Co. Ltd](#) 2019-nCoV Antibody Test (Colloidal Gold) (China NMPA EUA - Australia TGA - Brazil ANVISA - Singapore HSA - CE-IVD) [Contact](#)
- [InTec Products, Inc.](#) Rapid SARS-CoV-2 Antibody Test (CE-IVD) [Contact 1](#) [Contact 2](#)
- [InTec Products, Inc.](#) Rapid SARS-CoV-2 Antibody (IgM/IgG) (CE-IVD) [Contact 1](#) [Contact 2](#)
- [Jetta Labs LLP](#) OZO Diamond SARS-CoV2 (COVID-19) IgG/IgM Test (Latex Method) (CE-IVD) [Contact](#)
- [Jetta Labs LLP](#) OZO India SARS-CoV2 (COVID-19) IgG/IgM Test (Colloidal Gold Method) (CE-IVD) [Contact](#)
- [Jiangsu Bioperfectus Technologies Co. Ltd](#) PerfectPOC Novel Corona Virus (SARS-CoV-2) IgM/IgG Rapid Test Kit (CE-IVD) [Contact](#)
- [Jiangsu Bioperfectus Technologies Co. Ltd](#) PerfectPOC Novel Corona Virus (SARS-CoV-2) Ag Rapid Test Kit (CE-IVD) [Contact](#)
- [Jiangsu Superbio Biomedical Technology \(Nanjing\) Co., Ltd](#) SARS-CoV-2 (COVID-19) IgM/IgG Antibody Fast Detection Kit (Colloidal Gold) (US FDA EUA - CE-IVD) [Contact](#)
- [JinHuan Medical Instrument Co., Ltd](#) (COVID-19) IgM/IgG Antibody Fast Detection Kit (Colloidal Gold) (CE-IVD) [Contact](#)
- [Joinstar Biomedical Technology Co., Ltd](#) SARS-CoV-2 IgM/IgG Antibody Test (Colloidal Gold) (CE-IVD) [Contact](#)
- [JOYSBIO \(Tianjin\) Biotechnology Co., Ltd](#) COVID-19 IgG/IgM Rapid Test Kit (Colloidal Gold) (CE-IVD) [Contact](#)
- [JOYSBIO \(Tianjin\) Biotechnology Co., Ltd](#) COVID-19 (SARS-CoV-2) Antigen Rapid Test Kit (Colloidal Gold) (CE-IVD) [Contact](#)
- [JOYSBIO \(Tianjin\) Biotechnology Co., Ltd](#) COVID-19 Neutralizing Antibody Test Kit (Lateral Flow Rapid Test) (CE-IVD) [Contact](#)
- [Kephra Diagnostics](#) KDX Rapid SARS-CoV-2 Antigen Test (In development) [Contact](#)
- [Kephra Diagnostics](#) KDX COVID-19 IgG/IgM Rapid Detection Test Kit (In development) [Contact](#)
- [Koch Biotechnology \(Beijing\) Co., Ltd](#) SARS-CoV-2 Antigen Lateral Flow Assay (MHRA UK) [Contact](#)
- [KRISHGEN BioSystems](#) Human Anti-SARS-CoV-2 (Covid-19) IgG/IgM Rapid Test (CE-IVD) [Contact](#)
- [KRISHGEN BioSystems](#) Human Anti-SARS-CoV-2 (Covid-19) IgM Rapid Test (RUO) [Contact](#)
- [L&H Biotech Limited](#) COVID-19 Antigen Rapid Test (In development) [Contact](#)
- [Labnovation Technologies Inc.](#) COVID-19 (SARS-CoV-2) IgM/IgG Antibody Test Kit (CE-IVD) [Contact 1](#) [Contact 2](#)
- [Labtest Diagnostica SA](#) Anti COVID-19 IgG/IgM Rapid Test (Brazil ANVISA) [Contact](#)
- [Leadgene Biomedical, Inc.](#) Leadgene® SARS/SARS-CoV-2 Antigen Rapid Test Kit (In development) [Contact](#)
- [Leadgene Biomedical, Inc.](#) Leadgene® SARS/SARS-CoV-2 IgG/IgM Rapid Test Kit (In development) [Contact](#)
- [Lifeassay Diagnostics Pty Ltd](#) Test-it COVID-19 IgM/IgG Lateral Flow Assay (In development) [Contact](#)
- [LifeSensors, Inc.](#) COVID-19 IgG ELISA Detection Kit (RUO) [Contact](#)
- [Liming Bio-Products Co., Ltd](#) COVID-19 IgG/IgM Combo Rapid Test Device (CE-IVD) [Contact](#)
- [LOMINA AG](#) Fast COVID19 IgM/IgG Antibody Detection Kit (Colloidal Gold) (CE-IVD) [Contact](#)
- [Luminostics, Inc.](#) CLIP-COVID19 (smartphone-read out high sensivity antigen detection test) (In development) [Contact](#)

Search Website

<https://www.finddx.org/covid-19/pipeline/>

The Emergency Use Listing

SARS-CoV-2 Rapid Antigen Tests: progress of the active applications in the emergency use listing assessment pipeline

Product name	Product code(s)	Manufacturer name	Dossier review	QMS Desk Assessment
ESPLINE SARS-CoV-2	231906	Fujirebio, Inc	R	
BIOEASY Diagnostic kit for SARS-CoV-2 Ag (Fluorescence Immunochromatographic Assay)	YRLF04401025, YRLF04401050 and YRLF04401100	Shenzhen Bioeasy Biotechnology Co., Ltd	awaiting submission	awaiting submission
LumiraDx SARS-CoV-2 Ag Test	L0160001nnxxx	LumiraDx UK Ltd	awaiting submission	awaiting submission
SARS-CoV-2 Rapid Antigen Test	9327592190	Roche Diagnostics GmbH		
SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)	G10313	JOYSBIO (Tianjin) Biotechnology CO., LTD		

Progress of the active applications in the emergency use listing assessment pipeline.

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