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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

APPLICATION

A modern device for removing freckles, moles, discolorations and minor skin lesions allows for effective and safe treatments in home conditions. Thanks to plasma technology (or electrocoagulation - depending on the model), the precise tip acts on the selected skin area, removing unwanted changes without damaging the surrounding tissue.

The device is intended for home use only and is not intended for medical or commercial use and is not a substitute for medical therapy. The device is intended for indoor use in dry rooms. It is not intended for medical or therapeutic use or for commercial use.

For safety and CE certification reasons, the product may not be rebuilt or modified in any way. If the product is used for purposes other than those previously described, the product may be damaged. Improper use may also cause hazards such as short circuits, fire, electric shock, etc.

ACTIVATION

- Note: Packaging material must be completely removed from the product.
- Make sure that the device components are properly assembled.
- Before first use, charge the device with a USB cable. Connect the USB cable to a USB power source and connect the USB cable plug to the USB input port of the product.
- Select the appropriate tip (needle) depending on the type of treatment and gently screw it into the tip of the device.
- Turn on your device by holding the power button for 2-3 seconds until the LED screen lights up.
- Use the buttons to set your preferred power level – it is recommended to start with the lowest.
- Test the device on a small, less visible area of skin.
- Holding the device perpendicular to the skin, touch the tip to the selected lesions – do not press too hard.
- Once the treatment is complete, turn the device off by holding the power button and then gently unscrew the needle.
- Clean the tip according to the instructions and store the device in a dry place.

TIP

- The device must be checked for completeness of delivery and any visible damage.
- In the event of an incomplete delivery or if you notice damage due to faulty packaging or transport, please contact the service hotline.

TECHNICAL DATA

- Clear LCD display
- Power regulation
- Thin and thick needles included
- Charging via micro USB cable
- Power supply: battery
- Power: 1.9W

SAFETY INSTRUCTIONS

- This product may be used by children aged 8 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- Never immerse the device in water.
- For cleaning, use a damp cloth or a mild detergent.
- Do not use a damaged device.
- The product is not intended for use by children.
- Make sure that packaging materials are not left unattended. Children may start playing with them, which is dangerous.
- Protect the product from extreme temperatures, direct sunlight, strong vibrations, high humidity, moisture, flammable gases, vapors and solvents.
- Do not expose the product to mechanical stress.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - does not function properly, - has been stored for an extended period under unfavourable conditions or - has been excessively stressed during transport.

- It is forbidden to use the product if any of its parts are damaged. In the event of damage to the cable, it is forbidden to make repairs yourself.
- Do not disassemble the device yourself.
- Protect the product from moisture.
- Make sure not to damage the USB cable with sharp edges or hot objects. Fully unwind the USB cable before use.
- The included USB cable is only suitable for this product.
- Immediately turn off the product and remove the cable from the product if it smells burning or smoke is visible. Have the product checked by a qualified technician before using it again.
- If a USB power supply is used, the socket must always be easily accessible so that the USB power supply can be easily removed from the socket in an emergency. Please also observe the manufacturer's operating instructions.
- Before each charging process, check the USB cable and device for any damage. Never use or charge the device if you find any damage.
- Before use, make sure that the existing mains voltage matches the required operating voltage of the device.
- Use only as intended – for external use on skin only.
- Do not use on irritated, open, infected or wounded skin.
- Do not use on eyelids, around the eyes or on mucous membranes.
- Before first use, read the operating instructions and manufacturer's recommendations carefully.
- Always disinfect needles before and after use – use disposable or thoroughly sterilized tips.
- Do not use the device near water or with wet hands.
- Keep the device away from children and animals – it is not a toy.
- Do not perform the procedure in one place for too long to avoid burning the skin.
- Do not use on moles without prior consultation with a dermatologist.
- After the treatment, avoid sun exposure, do not use cosmetics with alcohol and do not irritate the treatment area.
- It is recommended to test on a small patch of skin before first use.
- Do not disassemble or modify the device.
- In case of redness, burning or allergic reaction, discontinue use and consult your doctor.
- Always use the product as intended.

BATTERIES/ACCUMULATORS SAFETY INSTRUCTIONS

- Batteries/rechargeable batteries should be kept out of reach of children. If swallowed, consult a doctor immediately!
- Disposable batteries must not be recharged. Batteries/accumulators must not be short-circuited and/or opened. This may cause overheating, fire or explosion.
- Never throw batteries/accumulators into fire or water.
- Batteries/accumulators should never be subjected to mechanical stress.
- Risk of acid leakage from batteries/accumulators.
- Avoid extreme conditions and temperatures that may affect batteries / accumulators, e.g. radiators / direct sunlight.
- If batteries/accumulators leak, avoid contact of skin, eyes and mucous membranes with chemicals! Immediately rinse affected areas with clean water and seek medical attention!
- Leaking or damaged batteries/rechargeable batteries can cause chemical burns when in contact with skin. Therefore, suitable protective gloves should be worn in such cases.
- The product has a built-in battery that cannot be replaced by the user.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT

Used electrical and electronic products, due to environmental protection, must not be thrown into household waste with regular municipal waste, but disposed of properly. Information on collection points and their opening hours is provided by the relevant office.

This product complies with the requirements of the relevant European and national directives that apply to it. The product meets European and national requirements on the safety of devices and products.

This product complies with the requirements of the relevant European and national directives that apply to it, RoHS.

Improper disposal of batteries/accumulators poses a risk to the environment!

Batteries/accumulators should not be disposed of with household waste. They may contain harmful heavy metals and should be treated as special waste. The chemical symbols for heavy metals are as follows: Cd = cadmium, Hg = mercury, Pb = lead.

Therefore, used batteries / accumulators should be taken to municipal hazardous waste collection points.

The built-in battery must not be disassembled for disposal. The product must be handed over in its entirety to a collection point for used electronics.

We reserve the right to make changes to the text, design and technical data of the product without notice.

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