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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

APPLICATION AND DESCRIPTION OF THE DEVICE

The electric foot massager with EMS stimulation function is an ideal solution for those looking for relaxation and relief after a long day. The device uses electrical impulses that stimulate muscles and improve blood circulation, providing a feeling of relaxation and regeneration.

The device is intended for home use only and is not intended for medical or commercial use and is not a substitute for medical therapy. The device is intended for use inside dry rooms. It is not intended for medical or therapeutic use or for commercial use.

For safety and CE certification reasons, the product may not be rebuilt or modified in any way. If the product is used for purposes other than those previously described, the product may be damaged. Improper use may also cause hazards such as short circuits, fire, electric shock, etc.

WARNING

Danger of improper use!

Using the device for purposes other than its intended purpose and/or using the device in any other way may involve various risks.

The device must only be used for its intended purpose.

The procedures described in this instruction manual must be followed.

Warning! The device should not be used by people with a pacemaker, insulin pump, malignant tumors, pregnant women or menstruating women.

FUNCTION BUTTONS

- "M" - button for changing the settings of individual massager functions
- "+" and "-" - increases or decreases the selected option
- "SSS" - turns the heating on or off, each time you press this button the heating level changes
- "I/O" - power on/off

COMMISSIONING/INSTALLATION INSTRUCTIONS

- Note: Packaging material must be completely removed from the product.
- Make sure the components are assembled correctly.
- Place the massager on a flat, stable surface – preferably the floor. Make sure there is a power outlet within reach.
- Connect the power supply cable to the massager and then the plug to the socket. (This step can be skipped if you have previously charged the device)
- Sit comfortably and place your feet on the massage panels (bare feet for better conduction of EMS impulses).
- Turn on the massager using the button on the touch panel or use the remote control.
- Select one of the available massage modes ("M" button)
- Adjust the pulse level using the "+" and "-" buttons.
- After use, turn off the device and unplug the power supply from the socket.
- Wipe the massager panels with a dry or slightly damp cloth

TIP

- The device must be checked for completeness of delivery and any visible damage.
- In the event of an incomplete delivery or if you notice damage due to faulty packaging or transport, please contact the service hotline.

TECHNICAL DATA

- Material of construction: plastic
- Remote control
- Adjusting the massager's inclination
- Stable base
- Wide range of massage options
- Battery or mains power supply 100-240V AC, 50/60Hz

SAFETY INSTRUCTIONS

- This product may be used by children aged 8 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- Children are not allowed to play with the product.

- Cleaning and maintenance may not be performed by children without supervision.
- Never immerse the device in water.
- For cleaning, use a damp cloth or a mild detergent.
- Do not use a damaged device.
- Make sure that packaging materials are not left unattended. Children may start playing with them, which is dangerous.
- Protect the product from extreme temperatures, direct sunlight, strong vibrations, high humidity, moisture, flammable gases, vapors and solvents.
- Do not expose the product to mechanical stress.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - does not function properly, - has been stored for an extended period under unfavourable conditions or - has been excessively stressed during transport.
- It is forbidden to use the product if any of its parts are damaged. In the event of damage to the cable, it is forbidden to make repairs yourself.
- Do not disassemble the device yourself.
- Protect the product from moisture.
- Immediately turn off the product and remove the cable from the product if it smells burning or smoke is visible. Have the product checked by a qualified technician before using it again.
- Before using the device, make sure that the mains voltage matches the required operating voltage of the device.
- The device should not be used by people with a pacemaker, insulin pump, malignant tumors, pregnant or menstruating women.
- Always use the product as intended.

ELECTRIC SHOCK HAZARD

- Connect the device only to a properly installed and grounded power outlet. The mains voltage must match the data on the nameplate of the device / mains adapter.
- Make sure that the power cable does not get wet or damp during operation. The cable must be routed so that it is not pinched or damaged.
- Keep the power cable away from hot surfaces.
- Do not perform any repairs to the device yourself.
- Any repairs may only be performed by service or by qualified specialists.
- Before cleaning the device, unplug the power supply from the wall outlet.
- Do not use the device in a damp or wet environment, such as a bathroom.
- Have a damaged plug or power cable repaired immediately by a qualified specialist or service center in order to avoid any danger.
- Do not insert any needles or sharp objects into the device.

BATTERIES/ACCUMULATORS SAFETY INSTRUCTIONS

- Batteries/rechargeable batteries should be kept out of reach of children. If swallowed, consult a doctor immediately!
- Disposable batteries must not be recharged. Batteries/accumulators must not be short-circuited and/or opened. This may cause overheating, fire or explosion.
- Never throw batteries/accumulators into fire or water.
- Batteries/accumulators should never be subjected to mechanical stress.
- Risk of acid leakage from batteries/accumulators.
- Avoid extreme conditions and temperatures that may affect batteries / accumulators, e.g. radiators / direct sunlight.
- If batteries/accumulators leak, avoid contact of skin, eyes and mucous membranes with chemicals! Immediately rinse affected areas with clean water and seek medical attention!
- Leaking or damaged batteries/rechargeable batteries can cause chemical burns when in contact with skin. Therefore, suitable protective gloves should be worn in such cases.

WARNING! INJURY RISK!

- This device must not be used to massage infants, small children, helpless persons or animals.
- Never use this device if you experience dislocation or pain in the area being massaged.
- Stop the massage immediately if it becomes uncomfortable or causes pain.
- The device has a hot surface during operation. People who are not sensitive to heat should exercise caution when using the device.
- Do not use the device on swollen, burned, irritated, diseased or injured areas of the skin or body.

- If you have any doubts regarding use, please consult your doctor.
- The device must not be used on the face (e.g. eyes), larynx, feet or other sensitive parts of the body.
- Do not cover the device (pillow, blanket, etc.) during use. Do not use the device near gasoline or other flammable materials.

WARNING – MATERIAL DAMAGE!

Use the appliance indoors only.

Only use the device with the supplied mains adapter.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT

Used electrical and electronic products, due to environmental protection, must not be thrown into household waste with regular municipal waste, but disposed of properly. Information on collection points and their opening hours is provided by the relevant office.

This product complies with the requirements of the relevant European and national directives that apply to it. The product meets European and national requirements on the safety of devices and products.

This product complies with the requirements of the relevant European and national directives that apply to it, RoHS.

Electric shock warning! Danger to life!

Improper disposal of batteries/accumulators poses a risk to the environment!

Batteries/accumulators should not be disposed of with household waste. They may contain harmful heavy metals and should be treated as special waste. The chemical symbols for heavy metals are as follows: Cd = cadmium, Hg = mercury, Pb = lead.

Therefore, used batteries / accumulators should be taken to municipal hazardous waste collection points.

The built-in battery must not be disassembled for disposal. The product must be handed over in its entirety to a collection point for used electronics.

We reserve the right to make changes to the text, design and technical data of the product without notice.

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